

Step One: Expose Fish

Expose Fish and Harvest Tissues




Step Two: RNA isolation

Isolate RNA with preferred method




Step Three: cDNA Synthesis

Synthesize cDNA using reverse transcriptase, Bi-NVT primer and Transcript Switching Oligo (TSO)



Step Four: Normalize

Denature and Reanneal, then bind to Streptavidin (SA)



Step Five: Fragment to 750bp

Cut with Fragmentase to 750bp




Step Six: Adaptor Ligation

End Repair, A-Tail and Ligate PCR Adaptors



Step Six: Final PCR to enrich and add biotin


Can be repeated indefinitely for as much probe as needed. Uses 5' Biotinylated primer necessary for DNA capture



Final Probe Product


Step One: Isolate gDNA

Isolate genomic DNA with magnetic bead protocol




Step Two: Tagmentation

Tagment to 1-4kb using R1 and R2 with Tn5 variant




Step Three: Elongate and add i7 indexes by PCR

PCR with P7-i7-R2 and R1 primers and Pool at the plate level




Step Four: Hybridize and Capture

Mix with biotinylated cDNA probes and capture with SA and elute ssDNA



Step Five: PCR to elongate and make dsDNA


Elongate with P7 and R1 primers to make ds



Step Six: Tagmentat to ~500bp


Fragment DNA with Tn5 to 500 bp for sequencing

Product:



Step Eight: Elongate and add i5 indexes by PCR

PCR with R1-i5-P5 primers and P7 primer.



Step 9: Sequence with Illumina

